

How To Measure


All forms/manikins in this catalog are copy written. Any copying of these forms/manikins is prohibited and violators will be prosecuted. Pivot Pedestal manikins are patented (6,458,434).

Before altering or prepping forms, please check hide's fit on form. Once a form has been altered or prepped there will be NO EXCHANGE OR REFUND GIVEN.

All measurements are made WITHOUT SKIN.
All measurements are taken with a 1/4" metal measuring tape.
Using other tape measuring devices may give you measurements that may not match the measurements of Classics Manikins.

Joe Coombs Classics sizing method for manikins may not be the same as other supply companies.

Primary measurements should be Size A and Size B.



Size A Tip of nose to corner of eye as illustrated.

Size B At beginning of swell from under chin as illustrated. Approximately 2" down from small of neck.

Size C At small of neck as illustrated

Turns are determined by the animal's right or left.

Use Size A and Size B measurements for Antelope and Caribou.

AUTO-SET™ ANTLERS- The quick and easy way to attach antlers to a Classics for a perfect fit.

A- Cut half-way between back of eye socket and base of burrs.

B- Cut 1" to 1 1/2" below back of skull.

Following these instructions will assure you a perfect fit on a CLASSICS.

AUTO-SET™ EYES- The eye sockets of these forms are preset at the proper angle forward and tilt downward, according to our study of live whitetails. The marks at the front and back of the eye socket allow you to line up the pupil with these marks and the eye is automatically set.

PLEASE NOTE: An Asterisk (*) will be found by the form number and noted at the bottom of a description for AUTO-SET™ eyes.

Measuring Life sized Mammals

E-N = Eyes to Nose

N-T = Nose to base of the tail

C = Circumference around the largest part of the mid section.

